


Auto Tilt Table Conble 104AD

COLOR



SPEC

Size(mm)	2000 x 700 x 480~930 (WxLxH)
Weight(kg)	120Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt



Features

- Low noise, low vibration electric motor
- Elevation controlled by electric motor(48~93cm)
- Tilt table angle adjusted by electric motor(0~90°)
- Foot plate angle adjusted by electric motor
- Built-in safety height sensor, Display current height
- No right/left moving in table at adjusting up/down
- with X type moving method
- One touch type moving & fixing device system
- Separable left and right Arm rests
- Digital Display, Angle & Time set, Automatic sequential, Alarm function