



2019
PRODUCTS
CATALOGUE





Electric training stairs

Smart Step

SPEC

Size(mm)	3100x1200x1500(WxLxH)
Main Power	AC 220, 60Hz
Start	Hand Switch
Components	Body, Hand Switch
Weight(kg)	170KG



Features

- ☐ Saving space by using a Parallels Bar and Stair gait trainer
- ☐ Gait train for critical patient and senior
- ☐ Low noise, low vibration electric motor
- ☐ Adjustable height(0~50cm)
- ☐ Simple control by hand switch

Gait Trainer

Walking Romo

SPEC

Size(mm)	1010x942x1635~2310(WxLxH)
Main Power	AC 220, 60Hz
Start	Hand Switch, Lifting Hand Switch, Tablet PC
Components	Body, Hand Switch, Lifting Hand Switch, Tablet PC, Harness
Weight(kg)	150KG



Features

- ☐ Low noise, low vibration electric motor
- ☐ Walking speed adjustment : 20-step
- ☐ Stride control : 50-step
- ☐ Exercise time set function
- ☐ Harness and lifting system
- ☐ Handswitch control programmed tablet PC
- ☐ Knee angle & ankle angle adjustment device
- ☐ Operate Active mod / Passive mode
- ☐ Individual management system

Leg Rehabilitation Equipment

Conble LS

COLOR



SPEC

Size(mm)	2100x660x580~825(WxDxH)
Weight(kg)	164Kg
Main Power	AC220, 60Hz
Start	Hand Switch
Drive	Electronic DC 24V
Components	Body, Hand Switch, Fixed belt



Features

- ☐ Sliding motion allowing for weight bearing exercise
- ☐ Slope standing rehabilitation
- ☐ The angle of slide table controlled by electric motor
- ☐ (0~43°)
- ☐ Elevation controlled by electric motor(58~82.5cm)
- ☐ Knee joint bending angle control function
- ☐ Separated foot plates' angle controlled by electric motor
- ☐ Slide mat bench stopper & protractor attachment
- ☐ Leg isolation device for CIMT

Auto Tilt Table


Conble 104AD

COLOR



SPEC

Size(mm)	2000 x 700 x 480~930 (WxLxH)
Weight(kg)	120Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor(48~93cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Foot plate angle adjusted by electric motor
- ☐ Built-in safety height sensor, Display current height
- ☐ No right/left moving in table at adjusting up/down
- ☐ with X type moving method
- ☐ One touch type moving & fixing device system
- ☐ Separable left and right Arm rests
- ☐ Digital Display, Angle & Time set, Automatic sequential, Alarm function

Auto Tilt Table



Conble 104D

COLOR



SPEC

Size(mm)	2000 x 700 x 480~930 (WxLxH)
Weight(kg)	120Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor(48~93cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Right/Left foot plates angle adjustment
 - manual fnc.
- ☐ Built-in safety height sensor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ One touch type moving & fixing device system
- ☐ Separable left and right Arm rests
- ☐ Digital Display, Angle & Time set, Automatic sequential, Alarm function

Auto Tilt Table


Conble 104A

COLOR



SPEC

Size(mm)	2000 x 700 x 480~930 (WxLxH)
Weight(kg)	120Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor(48~93cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
foot plates angle adjusted by electric motor
- ☐ Built-in safety height sensor
- ☐ No right/left moving in table at adjusting up/down with X
type moving method
- ☐ One touch type moving & fixing device system
- ☐ Separable left and right Arm rests

Auto Tilt Table


Conble 104

COLOR



SPEC

Size(mm)	2000 x 700 x 480~930 (WxLxH)
Weight(kg)	120Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor(48~93cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Right/Left foot plates angle adjustment – manual fnc.
- ☐ Built-in safety height sensor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ One touch type moving & fixing device system
- ☐ Separable left and right Arm rests

Auto Tilt Table

Conble 54F

COLOR



SPEC

Size(mm)	2000 x 700 x 540 (WxLxH)
Weight(kg)	65Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Low noise, low vibration electric motor
- ☐ Simple control by hand switch
- ☐ Fixed table height (54cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Set of 4 locking casters
- ☐ Right/Left foot plates angle adjustment – manual fnc.

Auto Tilt Table

Conble 54

COLOR



SPEC

Size(mm)	2000 x 700 x 540 (WxLxH)
Weight(kg)	60Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Low noise, low vibration electric motor
- ☐ Simple control by hand switch
- ☐ Fixed table height (54cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Set of 4 locking casters

Rehabilitation table for children

Conble CT

COLOR



SPEC

Size(mm)	2000 x 700 x 540 (WxLxH)
Weight(kg)	60Kg
Main Power	AC 220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Fixed belt(3EA)



Features

- ☐ Work therapy and tilt table for children
- ☐ Prone stander function
- ☐ Simple control by hand switch
- ☐ Fixed table height (55cm)
- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Right/Left foot plates angle adjustment
 - manual fnc.
- ☐ Therapy table attached(adjustable angle)
- ☐ Adjustable length for patient's height
- ☐ Crumpling prevention
 - Fixing ankle function Set of 4 locking casters

Electric Table

Conble NS1

COLOR



SPEC

Size(mm)	1950 x 700 x 490~930 (W x L x H)
Weight(kg)	80Kg
Main Power	AC220, 60Hz
Start	Hand Switch / 3sides bar Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, 3sides bar Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ Face hole for breathing
- ☐ One touch method moving & lock device attached
- ☐ Control with Hand switch & 3sides aluminum bar switch

Electric Table

Conble NS2

COLOR



SPEC

Size(mm)	1950 x 700 x 490~930 (W x L x H)
Weight(kg)	80Kg
Main Power	AC220, 60Hz
Start	Hand Switch / 3sides bar Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, 3sides bar Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ One touch method moving & lock device attached
- ☐ Control with Hand switch & 3sides aluminum bar switch
- ☐ Gas shock absorber adjusting angle in head part
- ☐ Face hole for breathing

Electric Table

Conble NS3

COLOR



SPEC

Size(mm)	1950 x 700 x 490~930 (W x L x H)
Weight(kg)	110Kg
Main Power	AC220, 60Hz
Start	Hand Switch / 3sides bar Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, 3sides bar Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor
- ☐ Gas shock absorber adjusting angle in head part
- ☐ Back and leg part angle adjusted by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ Face hole for breathing
- ☐ One touch method moving & lock device attached
- ☐ Control with Hand switch & 3sides aluminum bar switch


Electric Table Conble G1

COLOR



SPEC

Size(mm)	1950 x 700 x 460~900 (W x L x H)
Weight(kg)	68Kg
Main Power	AC220, 60Hz
Start	Foot Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Simple control by foot switch
- ☐ Elevation controlled by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ No scratch the floor by rail moving method
- ☐ Attached the wheel (2EA)
- ☐ Simple design & Reasonable price

Bobath, Voita Table

Conble NB1

COLOR



SPEC

Size(mm)	2000 x 1200 x 460~900 (W x L x H)
Weight(kg)	105Kg
Main Power	AC220, 60Hz
Start	Hand Switch / 3sides aluminum bar switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, 3sides bar Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ One touch method moving & lock device attached
- ☐ Control with Hand switch & 3sides aluminum bar switch

Bobath, Voita Table

Conble NB2

COLOR



SPEC

Size(mm)	2000 x 1200 x 460~900 (W x L x H)
Weight(kg)	110Kg
Main Power	AC220, 60Hz
Start	Hand Switch / 3sides aluminum bar switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, 3sides bar Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Elevation controlled by electric motor
- ☐ Head angle adjusted by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ One touch method moving & lock device attached
- ☐ Control with Hand switch & 3sides aluminum bar switch

Bobath, Voita Table

Conble GB1

COLOR



SPEC

Size(mm)	2000 x 1200 x 460~900 (W x L x H)
Weight(kg)	80Kg
Main Power	AC220, 60Hz
Start	Foot Switch
Drive	Electric motor DC 24V
Components	Body, Foot Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Simple control by foot switch
- ☐ Elevation controlled by electric motor
- ☐ No right/left moving in table at adjusting up/down with X type moving method
- ☐ No scratch the floor by rail moving method
- ☐ Attached the wheel(2EA)
- ☐ Simple design & Reasonable price

Shower Bed

Conble ST

SPEC

Size(mm)	Size(mm) : 2000 x 810 x 570~1020 (W x L
Weight(kg)	x H)
Main Power	Main Power : AC220, 60Hz
Start	Start : Hand Switch
Drive	Drive : Electric motor DC 24V
Components	Components : Body, Hand Switch



Features


- ☐ Electric shower bed for critical patient
- ☐ Low noise, low vibration electric motor
- ☐ Simple control by hand switch
- ☐ Elevation controlled by electric motor
- ☐ Open & Close rail handle
- ☐ Enabled manual adjusting angle of shower bed for drainage
- ☐ Quick drainage by one touch drain hole (2EA)
- ☐ Built-in rechargeable battery

Bed for patient / Tilt Table

Conble TB

SPEC

Size(mm)	2000 x 1000 x 540 (W x L x H)
Weight(kg)	120Kg
Main Power	AC220, 60Hz
Start	Hand Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Table, Fixed belt(3EA)



Features

- ☐ Tilt table angle adjusted by electric motor(0~90°)
- ☐ Head part angle adjusted by electric motor(0~80°)
- ☐ Standing up training on bed
- ☐ Attached side rail
- ☐ Attached the folding table
- ☐ Tilt table function

C-Arm Electric Table

Conble C-Arm

SPEC

Size(mm)	2000 x 600 x 700~950 (W x L x H)
Weight(kg)	65Kg
Main Power	AC220, 60Hz
Start	Hand Switch / Foot Switch
Drive	Electric motor DC 24V
Components	Body, Hand Switch, Foot Switch



Features

- ☐ Low noise, low vibration electric motor
- ☐ Foldable leg support attached
- ☐ Easy fixing & moving device attached
- ☐ 3D fluoroscopy (MAT : No X-ray transparent obstacle)
- ☐ Minimum height 70cm, Maximum height 95cm, adjustable
- ☐ C-arm Range of movement : 1400mm
- ☐ One touch method moving & lock device attached
- ☐ Tilting backward and forward

Unloading walker KU-200



SPEC

Size(mm)	2000 x 1050 x 1900~2200 (W x L x H)
Weight(kg)	40Kg
Start	Hand Switch
Components	Body, Hand Switch, Sling sheet

Features

- ☐ Simple control by hand switch
- ☐ Height control by electric motor
- ☐ Built-in rechargeable battery
- ☐ Available on a treadmill
- ☐ Lock feature wheels attached



Unloading walker KU-200W



SPEC

Size(mm)	1200 x 1050 x 1900~2200 (W x L x H)
Weight(kg)	40Kg
Start	Hand Switch
Components	Body, Hand Switch, Sling sheet

Features

- ☐ Measure a patient's body weight
- ☐ Simple control by hand switch
- ☐ Height control by electric motor
- ☐ Built-in rechargeable battery
- ☐ Available on a treadmill
- ☐ Lock feature wheels attached

Standing Table KS-500



SPEC

Size(mm)	950 x 800 x 1050~1250 (W x L x H)
Weight(kg)	70Kg
Start	Hand Switch
Components	Body, Hand Switch, Sling sheet

Features

- ☐ Work therapy and stand in rehabilitation table
- ☐ Simple control by hand switch
- ☐ Auto Height control
- ☐ Standing motor control function
- ☐ One touch method ; moving table and fixed device attached
- ☐ Birch wood table



Standing Table KW-600



SPEC

Size(mm)	970 x 900 x 600~1000 (W x L x H)
Weight(kg)	30Kg
Start	Hand Switch
Components	Body, Fixed belt, Knee support

Features

- ☐ Standing table and work table function
- ☐ Table angle adjustment by manual
- ☐ Table height adjustment by manual
- ☐ Spacious width provided for wheel chair user
- ☐ Knee support for standing
- ☐ Fixed type belt for standing
- ☐ Easy moving and fixing through wheels
- ☐ Birch wood table



92, Hakjang-ro 169beon-gil, Sasang-gu, Busan, Republic of Korea
E.meditec@meditec.co.kr T.051-304-1520 F.051-304-1521
Copyright(c) 2019 Kwang Won Meditec All Reserved.